# TEG<sup>®</sup> 6s

Hemostasis System

## QC Level 1

**REF** 07-650

**LOT** 0830-1301


**HAEMONETICS<sup>®</sup>**

CE

IVD

**Rx Only**

**CONTROL**

HAEMONETICS S.A.  
Signy Centre  
Rue des Fléchères 6  
1274 Signy-Centre  
Switzerland

© 2015, 2016 Haemonetics Corporation.

Haemonetics, TEG, RapidTEG, and Thrombelastograph  
are trademarks or registered trademarks of Haemonetics Corporation in the  
United States and/or other countries.

PlateletMapping is a registered trademark of Cora Healthcare, Inc.  
US Patent Nos. 7,261,861, 7,879,615, 8,236,568, 9,068,966

## **EN - English**

### **Proprietary Name and Established Name**

**Proprietary Name:** TEG® 6s System and the following cartridge and Quality Control (QC) material:

- CaCl<sub>2</sub> (Calcium Chloride) Cartridge and Level 1 QC material

**Established Name:** Whole Blood Hemostasis System

### **Intended Use**

The TEG® 6s Hemostasis Analyzer is a non-invasive diagnostic instrument designed to monitor and analyze the coagulation state of a blood sample in order to assist in the assessment of patient clinical hemostasis conditions. The TEG 6s analyzer is indicated for use with adult patients where an evaluation of blood coagulation properties is desired.

Results from the TEG 6s analyzer should not be the sole basis for a patient diagnosis, but should be evaluated together with the patient's medical history, the clinical picture and, if necessary, other coagulation tests.

The TEG 6s System is used to test the hemostasis properties of blood samples. The TEG analyzer performs internal QC checks during a pretest when the cartridge is inserted. This verifies that all electromechanical and pneumatic functions of the analyzer-cartridge combination are operating satisfactorily. Laboratories should perform QC testing as required in accordance with current CLIA Regulations regarding levels and frequency of control materials. Haemonetics provides QC Level 1 and Level 2 materials and Calcium Chloride cartridges to be used for this purpose.

### **Principles of the Procedure**

The TEG Hemostasis System whole blood coagulation control Level 1 contains a preparation of citrated animal plasma. It has been formulated specifically as a quality control for the TEG analyzer. This control contains no human material. *For in vitro diagnostic use. For professional use only.*

### **Reagents**

CaCl<sub>2</sub> (Calcium Chloride) and Level 1 QC Material

- Dried Calcium Chloride (adequate to reverse citration in the QC material).
- QC Level 1 Material, consisting of one 1 mL vial of lyophilized citrated bovine plasma and one 1 mL vial of diluent water for reconstituting the lyophilized material.

### **Warnings and Precautions**

#### **For In-Vitro Diagnostic Use**

#### **Hazards**

The QC Level 1 material is potentially hazardous in case of ingestion, inhalation, skin contact and eye contact. It contains animal source material; therefore, it should be treated as potentially infectious. Use universal precautions for handling potentially bio-hazardous material. It may cause adverse effects for the environment if discarded improperly.

## Sources of Reagent Error

- Sealed cartridges should be refrigerated (2° C to 8° C) in sealed foil pouches when not in use.

## Sources of Procedural Error

- Carefully follow the instructions for the preparation of the QC Level 1 control material. For correct reconstitution, the full vial of diluent water must be poured into the vial of lyophilized control. Failure to do so may result in inaccurate results.
- Timing - Do not use controls if they have been reconstituted for more than 2 hours

## Sample Preparation

The entire contents of the diluent water vial must be placed into the Level 1 vial for proper concentration, and thoroughly mixed.

## Storage Instructions

Calcium chloride cartridges and QC materials must be protected from temperature extremes above 35°C. Haemonetics recommends that the cartridges and QC materials be refrigerated at 2°C to 8°C after receipt and prior to use; however, storage of up to one month at temperatures not exceeding 25°C will not have an adverse effect on the labeled shelf life or expiry date of the product. Cartridges are for single use only and must be discarded after use.

## Analyzer

Refer to the *TEG® 6s User Manual*.

## Specimen Collection and Preparation for Analysis

No specimen collection is required for running this quality control test. No additives or preservatives are necessary to maintain the integrity of the material, but reconstituted QC Level 1 material must be used within 2 hours of reconstitution.

## Materials Provided

- Ten (10) Calcium Chloride cartridges, ten (10) QC Level 1 vials and ten (10) Water vials, each sealed with a desiccant pack in a foil pouch. Each cartridge contains four (4) channels of dried calcium chloride.
- Twelve (12) transfer pipettes

## Materials Required (but not supplied)

TEG® 6s Analyzer (Source: Haemonetics)

## Test Procedure

1. Remove one cartridge pouch from refrigerated storage.
2. Tear open the pouch, remove the cartridge and vials, and discard the pouch and desiccant pack.
3. Prepare the QC sample according to the following instructions:
  - a. Allow the control vial and the diluent water vial to sit for approximately 10 minutes to equilibrate to room temperature.
  - b. Make sure the lyophilized material is on the bottom of the control vial. You may need to tap the vial a few times.
  - c. Remove the seal and stopper of the control vial, taking care to avoid sharp metal edges.

- d. Slowly pour the contents of the diluent water vial into the control vial. Make sure no water drips out.
  - e. Re-insert the stopper into the control vial.
  - f. Hold the stopper in place, vigorously shake the control vial until fully reconstituted, and then let it stand for 5 minutes at room temperature.
  - g. Shake the control vial vigorously again and let stand 5 more minutes.
4. From the *Home* screen on the TEG analyzer, touch **new qc**.
  5. Insert the cartridge into the slot, as indicated on the screen, with the bar code on the left side.
  6. After the cartridge pretest has completed and you have verified that the assay is what you intended to run, touch **next**. If the cartridge is not the desired assay, touch **stop** and start over.
  7. On the *Test Information* screen, enter information for the test (if desired), and then touch **next**.
  8. Pipette the prepared QC sample into the cartridge sample port, filling up to or above the line marked on the cartridge.
  9. Touch **next**.  
The TEG analyzer starts the test, and the results are displayed as they become available.
  10. When the analyzer displays the "Remove cartridge" prompt, remove the used cartridge from the slot and immediately dispose of it.

## Results

No calculations are necessary to obtain test results; these are displayed on the TEG analyzer display, and may be printed if a printer is attached.

## Limitations of the Procedure

### Use of Results

- Results outside of the normal ranges listed below may indicate a problem with analyzer performance, or it may be caused by improperly prepared control material or product deterioration. It is recommended that any failed QC control test should be repeated. Repeated failures should be reported to the supervisor.
- Do not report any patient results until the control parameters are within their appropriate ranges.

## Expected Values

Expected values for test results are within the reference ranges that were established based on testing, following Clinical and Laboratory Standards Institute (CLSI) publication C28-A3c, *Defining, Establishing and Verifying Reference Intervals*. Each laboratory may establish its own level of performance to monitor quality assurance. Follow your institution's policy to establish mean target values and ranges for the quality controls.

## Reference Ranges for QC Level 1

Citrated Blood Parameter	Range
R (min)	0.7 – 1.7
K (min)	0.4 – 0.6
Angle (deg)	77 – 86
MA (mm)	54 – 67

## Specific Performance Characteristics

### Precision

The precision of the TEG System QC Calcium Chloride cartridge with Level 1 was evaluated according to CLSI EP5-A2. Testing was performed for within-device (20-day tests using QC Level 1), and multi-site reproducibility tests were performed across three sites, three lots, three operators and three analyzers. Results included Coefficient of Variance (CV) values for all precision tests. The CV of the test results was less than 10% for all parameters.

### Legal Notice

Licensed for sale or use only in in-patient hospitals and any laboratories that are primary providers of such hospitals.

### Packaging Symbols

Symbol	Description
	Manufacturer (address for)
	Authorized representative in the European community (address for)
	CE mark
	Catalogue number
	Batch code
	Control
	In vitro diagnostic device
	Consult instructions for use
	Do not reuse
	Temperature limitation for storage
	(Applies to USA only) USA Federal Law restricts the sale, distribution or use of this device to, by or on the order of a physician.

Symbol	Description
	Contains sufficient for <n> tests
	Use by <date>
	Do not use if package is damaged

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